



USER MANUAL

A Visual Guide For Supporting Application Usage

Version 1.0

Table of Contents

Welcome Screen	3
Language	3
Grid Size	3
Main Screen	4
Sentence	4
Folders	4
Tiles	4
Edit Mode	5
Tile Arrangement	5
Column Locking	5
Tile Removal	5
Add/Edit Tile	6
Tile Name	6
Tile Image	6
Background Color	7
Custom Voice	7
Tile Behavior	8
Add Tile Set	9
Add Tiles	10
Add/Edit Folder	11
Folder Name	11
Folder Icon	11
Settings	12


Welcome Screen

Your preferred language.




This can be changed anytime in the language section within **Settings**.

Your preferred grid size.

The grid can be further customized anytime within **Edit Mode** or reset to a default preset within **Settings**.



Main Screen




-  Speak the current sentence.
-  Remove the last word.
-  Clear the entire sentence.

Sentence

Folders

Tiles



-  Lock / unlock setup icons.
-  Customize tiles, folders, layout.
-  Application settings.

Edit Mode

Touch or swipe
to the left for
additional options.

Tile and folder
management
options.

Draggable
tiles. Swap with
other tiles or
move into
another folder.

Order folders
by moving them
up and down.

Scroll zone
to allow for
vertical
scrolling.

Tile visibility
adjustment
toggle.

Lock or unlock
column content in
place.

Remove tiles
by dragging
them here.




Add/Edit Tile #1

Name of the current tile.

Custom image option, allowing for the use of an image stored on your device.

Remove image from current tile.

Available images for the given tile name.




Add/Edit Tile #2

Preset colors

available for setting the tile background color to. Based on the currently configured color scheme.

Custom voice

gives the option of recording a phrase for the tile, which then replaces the system voice.



Custom colors


allows for color palettes to be saved, managed and easily applied to tiles.

Tile Behavior

Tile option sets the tile to behave in its default manner.


Folder option sets the tile to behave as a folder. Upon selecting, a list of available target folders will appear below.

Manage behavior button is the entry point to this dialog.



Add Tile Set


Tile sets that are available to be added. They can be scrolled through in the list or found by using the search field below. They generally contain 40 relevant tiles each.



The location by which you wish to add the selected tile set contents to. Either to your currently selected folder or as a newly created folder.

Add Tiles

Color to use for all of the newly added tiles.




Tile names to be added as a batch set, allowing for quicker page setup. They will be ordered left to right, top to bottom. Names should be comma separated.

Add/Edit Folder


The icon to associate the folder with.

Remove the current folder along with all tiles inside.



The name to associate this folder with. This is more applicable on larger devices, where extra horizontal space allows for the name to be optionally displayed.

Settings



Info dialogs that contain a brief description for what every setting controls within the application.